# REUSABLE ENDOSCOPE CLEANING BRUSHES

(ENGLISH & FRENCH IFU)

CHANNEL CLEANING BRUSHES  
8107, 8115, 8160, 8170, 8180, 8182,  
8190, 8195, 8701, 8703, 8704


SUCTION WELL CLEANING BRUSH  
8880

Instructions for Use **2**

Mode d'emploi **6**

Choose PriMed for the Perfect Swipe!

## DEVICE DESCRIPTION – Channel Cleaning Brush

Suction Well Cleaning Brush

Figure 1

Rx ONLY

EC REP

PRIMED INSTRUMENTS INC.  
1080 Tristar Drive, Unit 14  
Mississauga, Ontario, L5T 1P1  
Canada  
PHONE: 905-565-0565  
FAX: 905-565-0566  
www.primedENDO.com

ADVENA LIMITED  
Tower Business Centre  
2<sup>nd</sup> Flr., Tower Street  
Swatara, BKR 4013  
Malta

Instructions for Use  
W-7.5.3-W4 - RCB IFU  
EN & FR  
Revised & Approved & Effective  
2020/01/13  
© 2020 PriMed Instruments Inc.

## INTENDED USE

These instruments are intended to be used as an aid in cleaning and disinfecting flexible endoscopes. Channel Cleaning Brushes are used to scrub the inside of the Endoscope Channel and loosen any debris present. Suction Well Cleaning Brush is used to scrub the inside of Suction Well of the Endoscope and loosen any debris present. Do not use these instruments for any purpose other than their intended use.

## SIGNAL SYMBOLS

<b>WARNING</b>	Indicates a potentially hazardous situation, which could result in death or serious injury.
<b>CAUTION</b>	Indicates a potentially hazardous situation, which may result in minor or moderate injury.
<b>NOTE</b>	Indicates additional helpful information.

## CONTRAINDICATIONS

Do not use the cleaning brushes if the channel(s) are known to be occluded as this may result in damage to the channel. Never use excessive force when introducing or withdrawing the brush. Use of excessive force may damage the brush or the equipment being cleaned.

Follow all cleaning instructions of the endoscope's Original Equipment Manufacturer. Do not use a specific cleaning brush in a specific channel if use is not recommended by the endoscope's Original Equipment Manufacturer.

## CAUTIONS

- Do not use Reusable Cleaning Brushes for cytological tissue sampling or any other diagnostic or therapeutic purposes. Patient injury, cross-contamination, equipment damage may occur.
- Never coil brushes into a diameter smaller than 20cm. The brushes can be damaged if coiled too tightly.
- Only trained and qualified personnel should be given the responsibility to reprocess endoscopes and their accessories.
- Ensure that the brush size is compatible with the Endoscope Channel & Suction Well.
- When using detergent and disinfection solution, always follow the manufacturer's recommendations regarding temperature, soaking time and concentration. Always use neutral pH solution for cleaning and disinfection. Subjecting brushes to highly acidic or alkaline solution for a prolonged period can lead to damages.
- Never use excessive force when introducing or withdrawing the brush. If a brush encounters significant resistance or blockage in the endoscope channel, do not force it through the endoscope. Use of excessive force may damage the endoscope channel.
- Do not rotate, twist or torque the brush while it is inside the channel.
- Do not use damaged brushes. See Figures 2-8. If a brush was damaged during preparation, inspection, use or reprocessing then it should be discarded to avoid damaging the endoscope.
- Perform cleaning and disinfection with both the brush and the endoscope completely submerged to prevent spattering.
- Never transfer a contaminated brush between endoscopes; cross-contamination may occur resulting in patient injury.
- Immediately after use, the brush must be cleaned to remove mucus, tissue and blood. Do not let a used brush dry before cleaning. Failure to properly clean the brush will compromise sterility and may cause the brush to fail. Inadequate cleaning of the brush will void the warranty, since accumulation of dried material on the brush will cause product deterioration.
- Cleaning used brushes with automated reprocessing system is NOT recommended.
- Sterilization must be performed via Steam Sterilization (Autoclave) or High Level Disinfection using only neutral pH solution. The use of automated cleaner/sterilizer system is NOT recommended. The use of Hydrogen Peroxide or Peracetic Acid as sterilizing agents is NOT recommended.

## WARNINGS

- Patient debris and reprocessing chemicals are hazardous. Wear protective equipment at all time.

## LIMITATIONS ON REPROCESSING & END OF DEVICE LIFE

Repeated processing has minimal effect on the cleaning brush. End of life is normally determined by wear and damage due to repeated use. Follow inspection instructions closely to identify brushes that no longer can fulfill their intended use. Failed inspection is an indication that used/reprocessed brush has reached its end of life. A brush, past its useful life, should be treated as hazardous waste and disposed of in accordance with hospital/clinic protocol.

## EQUIPMENT REQUIRED

- Preparation, Inspection and Use: Endoscope & reusable brush.
- Cleaning: Ultrasonic cleaner, container, enzymatic detergent, cleaning brush, sponge & lint-free cloth.
- Sterilization: Steam sterilizer (autoclave) or high level disinfectant.

## PREPARATION

- Before brushing, follow instruction manual provided by the endoscope manufacturer regarding endoscope preparation, inspection and cleaning.

## INSPECTION

- Remove the brush from the package and carefully unwind the main body from the proximal (handle) end.
- Inspect the brush between uses and replace when worn, frayed, bent, kinked, missing solder or otherwise damaged. See Figures 2-8. This inspection is the responsibility of the user. Worn out bristles are ineffective in cleaning, and damaged brushes may damage the endoscope.
- Do not repair damaged brushes. Use of repaired brushes may damage the endoscope.


Figure 2 Missing bristles (Gap)


Figure 3 Bent (Set) bristles


Figure 4 Loose (Wild) bristles


Figure 5 Permanent bristles adhesion

Figure 6 Damaged/kinked shaft


Figure 7 Brush sections that may detach


Figure 8 Brush head with missing solder

## USE

### CHANNEL CLEANING BRUSH

While the brush and the endoscope are submerged in a cleaning solution, advance the brush into the channel of the endoscope in short strokes to avoid kinking until the brush head emerges at the other end of the channel. With the brush still submerged, scrub off all debris on the bristles before reversing its direction in the channel. If using a double ended cleaning brush, make sure to scrub off all debris before pulling the other end of the brush through the channel.

### SUCTION WELL CLEANING BRUSH

While the brush and the endoscope are submerged in the cleaning solution, gently clean the well openings. Avoid twisting or torquing the brush handle.

## REPROCESSING

## CLEANING

- Prepare a large container of clean warm tap water and the enzymatic detergent (e.g. Enzol<sup>®</sup>) solution prepared according to manufacturer's instruction. Coil and place the brush into the solution, distal end first. Ensure the brush is completely immersed. Allow the brush to soak for a minimum of one (1) minute. Avoid prolonged soaking (e.g. overnight).
- Use a sponge or a lint-free cloth to thoroughly scrub the brush to remove adhering debris. Using finger tips, thoroughly scrub the brush bristles. If the brush is not thoroughly cleaned, residual organic debris may solidify and it may be difficult to properly disinfect the brush. Rinse the brush in a large container with ample amounts of clean warm tap water for a minimum of one minute.
- Proceed with the ultrasonic cleaning. Fill the ultrasonic tank with clean warm tap water and an enzymatic detergent (e.g. Enzol<sup>®</sup>) solution mixed according to manufacturer's instruction. Completely immerse the brush in the ultrasonic cleaner tank. Maintain the 20cm minimum coil diameter. Clean the brush in the ultrasonic cleaner with a minimum 40 kHz frequency output for five (5) minutes. If debris is still present, continue with longer ultrasonic cleaning. Rinse thoroughly under ample amounts of clean, running tap water for at least one (1) minute to remove any residual cleaning solution. Wipe dry with clean, dry, lint-free cloth.

## DISINFECTION

- Disinfect the brush in a High-Level Disinfectant (e.g. Cidex<sup>®</sup>) solution according to manufacturer's instructions regarding concentration, temperature and cycle time.
- While the brush is completely immersed, use a lint-free cloth to remove bubbles from outer surface of the brush. Use the gloved fingers tips to remove bubbles from the brush bristles.
- Rinse thoroughly with ample amounts of clean sterile water for at least one (1) minute. Repeat rinsing procedure TWO (2) additional times, for a total of THREE (3) rinses. Use fresh sterile water for each rinse.
- Use a lint-free cloth to thoroughly dry the brush.

NOTE: Contact PriMed Instruments for other brands of disinfectant that is compatible with our reusable cleaning brushes.

## STERILIZATION:

Thorough manual cleaning prior to sterilization is essential to achieve sterility. Organic debris on the brush and residual moisture will inhibit sterilization.

Place the reusable brush in a suitable container or package.

## CAUTIONS

- Do not crush bristle during sterilization.
- Follow the sterilizer manufacturer's instruction and do not exceed the parameters given in the table below.

Recommended autoclave cycles:

Sterilizer	Temperature	Exposure Time
Prevacuum	132°C (270°F)	4 minutes

NOTE: Contact PriMed Instruments for other acceptable parameters.

## STORAGE

Reusable Cleaning Brushes should be stored in a sterile package at room temperature in a dry clean area away from direct sunlight.

## PRODUCT DISPOSAL

Used cleaning brush may be a potential biohazard. Handle and dispose of in accordance with the hospital/clinic protocol and applicable regulations

## SYMBOLS MEANING

CE Mark.

Quantity of Units in Case.

Product Number.

Quantity of Units in Dispenser.

This product is not made with LATEX.

This product is not made with PHTHALATE.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Visit us online at [www.primedENDO.com](http://www.primedENDO.com) for full list of symbol glossary.

NOTE: Contact PriMed Instruments if you have any questions of comments regarding the content of this document.

If two years have elapsed between the issued date for these instructions and product use, the user should contact Primed Instruments to determine if additional information is available.

Enzol<sup>®</sup> is a registered trademark of Johnson & Johnson Inc.

Cidex<sup>®</sup> is a registered trademark of Johnson & Johnson Inc.

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

**ADVENA LIMITED**  
Tower Business Centre  
2<sup>nd</sup> Flr., Tower Street  
Swatar, BKR 4013  
Malte

Instructions for Use  
**W-7.5.3-W4 - RCB IFU**  
**EN & FR**  
Revised & Approved & Effective  
**2020/01/13**  
© 2020 PriMed Instruments Inc.

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