**PRIMED INSTRUMENTS INC.**

1080 Tristar Drive, Unit 14  
Mississauga, Ontario L5T 1P1, Canada

Phone: 905-565-0565

Fax: 905-565-0566

Email: [info@primedENDO.com](mailto:info@primedENDO.com)

TYPE: INSTRUCTIONS FOR USE

ID: W-7.5.3-W1A

**ADDENDUM to IFU W-7.5.3-W1**  
**for MultiCROC® (aka CROC® Multiple Sample Biopsy Forceps) models.**

**IMPORTANT:** READ AND FOLLOW ALL THE INSTRUCTIONS IN THE IFU BOOKLET W-7.5.3-W1 AND THE INSTRUCTIONS IN THIS ADDENDUM BEFORE USING THE MULTICROC® FORCEPS.

**INTENDED USE**

The CROC® Multiple Sample Biopsy Forceps (only model numbers starting with “MC”) have the capability to collect multiple biopsies in a single pass of the device through the channel of the endoscope (up to 5 tissue samples per single pass). Total number of biopsies collected will depend on tissue type and user’s preference.

**CAUTIONS**

- Do not attempt to pass the instrument when the scope is excessively articulated. Doing so will damage the endoscope
- MultiCROC® biopsy forceps jaws are designed to facilitate collection and storage of multiple biopsies. During each subsequent biopsy, tissues samples are pushed toward the back to fill-up the jaw’s elongated cups. Always monitor amount of samples stored inside the jaw’s cups and do not overload them. Overloading may result in crush artifacts and lost samples.

**WARNING**

- MultiCROC® biopsy forceps jaws have a smaller aperture, when opened, and a longer distal end than standard biopsy forceps. Do not apply excessive force to advance the forceps into the tissue. Applying excessive force may result in perforation of the biopsy site.
- The jaws are designed to travel about 1-2mm back during closing of the device. If necessary, compensate for this movement by gently forwarding the biopsy forceps toward the biopsy site. Make sure to apply a sufficient force to close the jaws completely. Maintain pressure to keep the jaws closed to retain the sample(s) until the next biopsy is taken.

**USE**

- To obtain subsequent samples, it is not necessary to withdraw the forceps from the endoscope. Up to five biopsy samples may be taken without withdrawing the forceps from the scope.  
**NOTE:** When performing biopsy in confined spaces such as in the Esophagus, pull jaws closer to tip of endoscope to avoid potential obstruction to jaws during opening. The Jaws may not open fully unless a slight continuous pressure is applied to the handle prior to taking the sample.  
**NOTE:** Do not drag open jaws loaded with harvested samples against the surrounding organ cavity and do not open or close the jaws rapidly. Doing so may result in a loss of tissue samples.  
**NOTE:** If a device is to be reused on the same patient during the same procedure, be sure that all biopsy samples are removed prior to reuse. Failure to do so may result in a misdiagnosis.

# **CROC<sup>®</sup>**

## **SINGLE-USE NON-THERMAL BIOPSY FORCEPS**

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CE 1639

## INTENDED USE

PriMed biopsy forceps are used to sample lesion(s) of the gastrointestinal track, such as polyps and suspected malignancies. CROC Biopsy Forceps are a single use device and should be discarded after each procedure. The physician is responsible for procedural preference and sample preparation, since these usually vary greatly. Do not use this device for any purpose other than the stated intended use.

## SIGNAL SYMBOLS

<b>WARNING</b>	Indicates a potentially hazardous situation, which could result in death or serious injury.
<b>CAUTION</b>	Indicates a potentially hazardous situation, which may result in minor or moderate injury.
<b>NOTE</b>	Indicates additional helpful information.

## POTENTIAL COMPLICATIONS

Potential complications associated with gastrointestinal endoscopy include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

## CONTRAINDICATIONS

Contraindications include those specific to any endoscopic procedure, as well as those specific to biopsy.

## CAUTIONS

Do not coil biopsy forceps into a diameter smaller than 20cm.  
Do not excessively twist (torque) the device during use or handling.  
Do not use excessive pressure when actuating the forceps. Excessive pressure will cause damage to the instruments.  
Keep forceps jaws closed during insertion into, traversal through, and removal from the endoscope channel.  
Forceps are compatible with the minimum channel size printed on label. Using forceps with incompatible channel size will damage endoscope.  
Special Caution should be taken when handling instruments that have a needle.

## REUSE WARNINGS

Instrument is sterilized with Ethylene Oxide (EO). Do not use if sterile package is punctured. If damage is found call PriMed Instruments Inc. or the authorized distributor for a return authorization.

Biopsy forceps is for SINGLE USE ONLY. Do not reprocess, re-sterilize or reuse. Reprocessing, re-sterilization or reuse may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reprocessing, re-sterilization or reuse may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

## WARNINGS

Biopsy forceps should only be used by, or under the supervision of, physicians trained in endoscopic tissue biopsy.

The CROC biopsy forceps should only be used to biopsy tissue where possible hemorrhage will not present a danger for patients. Adequate plans for management of potential bleeding and appropriate airway management should be in place.

Please follow inspection instructions below before using instruments. Failure to follow instructions may lead to serious medical consequences.

## STORAGE

Biopsy forceps should be stored in the sterile package at room temperature in a clean and dry area.

## EQUIPMENT REQUIRED

1. Endoscope.
2. Fixative & biopsy sample container.

## INSPECTION

1. Verify the expiration date on the package label. Do not use expired instrument.
2. Remove the forceps from the package and carefully unwind the sheath from the proximal (handle) end. Remove the protective guard from the jaws.
3. Hold the distal (jaw) end of the forceps in one hand and gently operate the jaws by working the finger spool. Normally, forceps should open and close easily with two loops of approximately 20cm in diameter. **CAUTION:** Operating the forceps while the unit is still tightly coiled may lead to damage.
4. Inspect for any functional abnormalities and any kinks or bends in the outer sheath. **CAUTION:** If any irregularities are noted, DO NOT USE the instruments. Call PriMed Instruments Inc. or the authorized distributor for a return authorization. Do not repair instruments. Attempts to repair nonfunctional or damaged instruments may induce dangerous failure modes during use.

## USE

1. First, verify that the channel size is adequate to accommodate the instrument with the jaws closed. The forceps jaws should be closed when introduced into the endoscope channel. Advance the instrument slowly down the channel, using small strokes, until the jaws can be visualized at the end of the scope. **NOTE:** It may be necessary to reduce the angle on the angulated portion of the scope to allow for ease of passage of the forceps. Do not force forceps through the endoscope channel.
2. To obtain a biopsy specimen, gently open the jaws and advance the forceps to the biopsy site. If needed, articulate the end of the endoscope to direct the forceps to the biopsy site. Close the jaws completely to cut cleanly through the tissue and gently pull the forceps back. **NOTE:** Do not apply excessive pressure when closing forceps jaws. Applying excessive force to the handle spool or thumb ring may damage the instrument.
3. Keep the jaws closed to retain the sample. Slowly withdraw the forceps. **NOTE:** If the biopsy sample is too large to close the jaws sufficiently, withdraw the forceps as far as the tip of the endoscope, and withdraw the scope and biopsy forceps at the same time.

## PRODUCT DISPOSAL

**WARNING:** These products may be a potential biohazard. Treat used forceps as hazardous waste and dispose of in accordance with hospital/clinic protocol.

## SYMBOLS MEANING

Quantity of Units in Case.

Quantity of Units in Dispenser.

CE Mark.

Product Number.

This product is not made with Latex.

This product is not made with Phthalate.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Visit us online at [www.primedENDO.com](http://www.primedENDO.com) for full list of symbol glossary.

NOTE: Contact PriMed Instruments if you have any questions or comments regarding the content of this document.

If two years has elapsed between the issued date for these instructions and product use, the user should contact Primed Instruments to determine if additional information is available.

CROC<sup>®</sup> is a registered trademark of PriMed Instruments Inc.

This product is protected by one or more of the following Patents;  
U.S. Patent No. 6,309,404; 7,171,839  
Canadian Patent No. 2,094,529  
International Patents  
Other Patents pending

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© 2020 PriMed Instruments Inc