

With **SmartSTOP**[™] - engineered for prevention and safety.


✔ Prevents Accidental Jamming

SmartSTOP[™] is a *Calibrated, Color-coded safeguard* that prevents accidental insertion into undersized channels during reprocessing. It protects endoscopes from costly damage while streamlining SPD/MDR technicians' workflows.

Yellow SmartSTOP[™]
blocks entry into
channels smaller
than 2.8 mm



Blue SmartSTOP[™]
blocks entry into
channels smaller
than 2.0 mm





With **SmartSTOP**[™] - engineered for prevention and safety.

✓ Instantly **Detects**  **Channel Flaws & Obstructions**

SmartSTOP[™] delivers tactile feedback during brushing, alerting SPD/MDR technicians to hidden channel flaws. Endoscopes with flaws can be quickly identified and serviced before use, safeguarding patients, saving time, resources and eliminating repeated procedures.

- Kinked or collapsed channels



- Obstructions from retained accessories (e.g., hemoclips)



primedENDO.com

TO ORDER CONTACT

 sales@primedENDO.com

 Toll Free: **877-565-0565**



Patents Pending. Manufactured in Canada under ISO13485, CE Mark, Health Canada & FDA listed.

Declaration of Performance and Validation data available upon request.

Printed in Canada LT092025