



NoSpill™

Prevents Formalin Spills

Eliminate:

- ⊘ Dangerous formalin spills
- ⊘ Lost tissue specimens
- ⊘ Tedious clean-ups



Various hole sizes available for different containers.



primedENDO.com

SIZING CHART

NS30204-100

NS45204-100

NS55204-100

Match your container
to check which size
is right for you.



NoSpill™

Prevents Formalin Spills

The International Agency for Research on Cancer classifies formaldehyde as a human carcinogen.

Formalin spills are dangerous and have to be reported to Occupational Health Departments.

NoSpill™ is a simple and effective solution to minimize and prevent formalin spills in areas where tissue samples are harvested and processed, such as Endoscopy Units, ORs and Pathology Departments.

EASY TO USE!

1

UNFOLD (OPEN)
THE BOX FROM
THE FLAT STATE



2

PRESS IN TABS
ON BOTH SIDES



3

INSERT CONTAINER
INTO CUTOUT



PART NUMBER	HOLE DIAMETER	UNITS PER BOX
NS30204-100	30 mm	100 pcs
NS45204-100	45 mm	100 pcs
NS55204-100	55 mm	100 pcs

Read Instructions for Use Pamphlet before using.

ORDERING INFORMATION

Email: sales@primedENDO.com

Phone Toll Free: **877-565-0565**

Fax Toll Free: **877-565-0566**

Manufactured by:
PriMed Instruments Inc.
1080 Tristar Drive, Unit #14
Mississauga, ON L5T 1P1
Canada

905-565-0565

905-565-0566



primedENDO.com

Patents Pending / ©2019 PriMed Instruments Inc. Manufactured under ISO 13485 & MDSAP Quality Systems. CE Mark, Health Canada & FDA listed.