



CHOOSE PRIMED FOR THE PERFECT SWIPE!

REUSABLE CHANNEL CLEANING BRUSH
PN# 8180



REUSABLE SUCTION WELL CLEANING BRUSH
PN# 8880

See specifications

CHOOSE PRIMED FOR THE PERFECT SWIPE!

Consider some of these innovative features:

- Stainless steel bullets minimize the wear of the endoscope channels
- Ergonomic design reduces user's fatigue and improves productivity
- Very smooth stainless steel coil sheath
- Stainless steel cable inside brush eliminates coil sloppiness
- Autoclavable and validated for cleaning efficacy*, sterilization parameters** and multiple reprocessing

SPECIFICATIONS FOR REUSABLE ENDOSCOPE CLEANING BRUSHES

Reusable Channel Cleaning Brush

PRODUCT NUMBER	DESCRIPTION OR FOR USE IN:	BRUSH HEAD DIAMETER	CATHETER DIAMETER	WORKING LENGTH	HANDLE COLOR	PACKAGING OPTIONS
8107	Cystoscope	3-5mm	1.6mm	70cm	White	Pouch of 1
8115	Bronchoscope	3-5mm	1.6mm	115cm	White	Pouch of 1
8160	Slim Gastroscope/ Ped	3-5mm	1.6mm	160cm	White	Pouch of 1
8170	Colonoscope/ Gastroscope	7.0- 9.0mm	1.8mm	240cm	White	Pouch of 1
8180	Colonoscope/ Gastroscope	3-6mm	1.8mm	240cm	White	Pouch of 1
8182	Colonoscope/ Gastroscope	3-6mm Double Ended	1.8mm	240cm	N/A	Pouch of 1
8190	Sigmoidoscope	3-6mm	1.8mm	120cm	White	Pouch of 1
8195	Colonoscope	3-6mm	1.8mm	250cm	White	Pouch of 1
8701	Wire Sheath, Micro Brush	2.3mm	0.5mm	95cm	White	Pouch of 1
8703	Wire Sheath, Micro Brush	2.3mm	0.5mm	165cm	White	Pouch of 1
8704	Wire Sheath, Micro Brush	2.3mm	0.5mm	120cm	White	Pouch of 1

Reusable Suction Well Cleaning Brush

PRODUCT NUMBER	DESCRIPTION OR FOR USE IN:	BRUSH HEAD DIAMETER	WORKING LENGTH	HANDLE COLOR	PACKAGING OPTIONS
8880	Suction Well	11mm	5.5cm	White	Pouch of 10 or 25

Other sizes and custom configurations available.

*AAMI TIR30:2011 post-processing TOC level <2.2µg/cm² and Protein level < 6.4µg/cm²

**Refer to IFU



ORDER TODAY!

Email: sales@primedENDO.com
Phone Toll Free: 877-565-0565
FAX Toll Free: 877-565-0566
Online: www.primedENDO.com



PriMed Instruments Inc.